


[KO Validated] Anti-Cdkn1B/P27Kip1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P46527

Item Number: CM0002905



Introduction

Validated knockout (KO) rabbit monoclonal antibody targeting CDKN1B/p27KIP1, suitable for Western blot, IHC-P, and ELISA in human, mouse, and rat samples. Ideal for cell cycle regulation studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] CDKN1B/p27KIP1 Rabbit mAb
Remarks/Alias	KIP1; MEN4; CDKN4; MEN1B; P27KIP1; CDKN1B/p27KIP1
Species	Human
GeneID (Human)	1027
GeneID	1027
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 99-198 of human CDKN1B/p27KIP1 (P46527).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P46527
Protein Weight	22kDa
Shipping	Ice bag