


[KO Validated] Anti-Cdkn2A/P16Ink4A Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P42771 / Q8N726

Item Number: CM0016859



Introduction

KO-validated rabbit polyclonal antibody targeting CDKN2A/p16INK4a, a tumor suppressor inhibiting CDK4/6. Suitable for WB, IF/ICC, and ELISA in human samples. Research-use IVD raw material. SWISS: P42771 / Q8N726.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] CDKN2A/p16INK4a Rabbit pAb
Remarks/Alias	ARF; MLM; P14; P16; P19; CMM2; INK4; MTS1; TP16; CDK4I; CDKN2; INK4A; MTS-1; P14ARF; P19ARF; P16INK4; P16INK4A; P16-INK4A; 4a
Species	Human
GeneID (Human)	1029
GeneID	1029
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-156 of human CDKN2A/p16INK4a (NP_000068.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P42771/Q8N726
Protein Weight	8kDa/11kDa/12kDa/13kDa/16kDa/17kDa
Shipping	Ice bag