

[KO Validated] Anti-Cetn2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P41208

Item Number: CM0026721



Introduction

Rabbit polyclonal antibody targeting CETN2/Centrin-2, validated for WB, IF/ICC, and ELISA. Reacts with human, mouse, and rat samples. Essential for centrosome duplication, DNA repair, and mRNA export studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] CETN2 Rabbit pAb
Remarks/Alias	CALT; CEN2; N2
Species	Human
Gene ID (Human)	1069
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-172 of human CETN2 (NP_004335.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P41208
Protein Weight	20kDa
Shipping	Ice bag