

[KO Validated] Anti-Chchd2 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Y6H1

Item Number: CM0019690



Introduction

[KO Validated] Rabbit polyclonal antibody against CHCHD2, a mitochondrial transcription factor implicated in Parkinson's disease. Validated for WB, IHC-P, IF/ICC, ELISA, with cross-reactivity to human, mouse, and rat.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] CHCHD2 Rabbit pAb
Remarks/Alias	MNRR1; NS2TP; MIX17B; PARK22; C7orf17; D2
Species	Human
GeneID (Human)	51142
GeneID	51142
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 75-145 of human CHCHD2 (NP_057223.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y6H1
Protein Weight	16kDa
Shipping	Ice bag