


[KO Validated] Anti-Chd1 Rabbit Polyclonal Antibody (Pab) For Wb, Ihc-P, Elisa - O14646

Item Number: CM0027705



Introduction

High-quality rabbit polyclonal antibody against human CHD1, validated for Western blot, IHC-P, and ELISA applications. Recognizes human, mouse, and rat CHD1. Suitable for chromatin remodeling, transcription, and splicing research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] CHD1 Rabbit pAb
Remarks/Alias	CHD-1; PILBOS; D1
Species	Human
GeneID (Human)	1105
GeneID	1105
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1600 to the C-terminus of human CHD1 (NP_001261.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14646
Protein Weight	197kDa
Shipping	Ice bag