


KO Validated Anti-Ciap1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q13490

Item Number: CM0003452



Introduction

KO validated rabbit monoclonal antibody targeting human cIAP1 (Q13490), suitable for WB, IHC-P, and ELISA. Cross-reacts with human and rat. Ideal for apoptosis, NF-kB signaling, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] cIAP1 Rabbit mAb
Remarks/Alias	API1; MIHB; HIAP2; RNF48; cIAP1; Hiap-2; c-IAP1; P1
Species	Human
GeneID (Human)	329
GeneID	329
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human cIAP1 (Q13490).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	Q13490
Protein Weight	70kDa
Shipping	Ice bag