


[KO Validated] Anti-Cirbp Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q14011

Item Number: CM0027336



Introduction

KO validated rabbit polyclonal antibody against cold-inducible RNA-binding protein (CIRBP). Suitable for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Stress response research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] CIRBP Rabbit pAb
Remarks/Alias	CIRP; BP
Species	Human
GeneID (Human)	1153
GeneID	1153
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-172 of human CIRBP (NP_001271.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14011
Protein Weight	19kDa
Shipping	Ice bag