

[KO Validated] Anti-Csde1 Rabbit Polyclonal Antibody For Wb, Elisa - O75534

Item Number: CM0022933



Introduction

Rabbit polyclonal antibody against human CSDE1 (UniProt O75534), validated for Western blot and ELISA applications. Reacts with human and mouse samples. Ideal for studying RNA-binding protein CSDE1 involved in mRNA turnover and stress granule formation.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] CSDE1 Rabbit pAb
Remarks/Alias	UNR; D1S155E; E1
Species	Human
GeneID (Human)	7812
GeneID	7812
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 369-600 of human CSDE1 (NP_009089.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O75534
Protein Weight	89kDa
Shipping	Ice bag