


[KO Validated] Anti-Csde1 Rabbit Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - O75534

Item Number: CM0030220



Introduction

KO-validated rabbit polyclonal antibody against CSDE1 (UNR), suitable for WB, IF/ICC, IP, ELISA. Detects human and mouse CSDE1. Ideal for translational regulation and stress granule research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] CSDE1 Rabbit pAb
Remarks/Alias	UNR; D1S155E; E1
Species	Human
GeneID (Human)	7812
GeneID	7812
Immunogen	Recombinant protein corresponding to amino acids 369-600 of human CSDE1 (NP_009089.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	O75534
Protein Weight	89kDa
Shipping	Ice bag