


KO Validated Anti-Cyclin B1 Rabbit Monoclonal Antibody For Wb And Elisa - P14635

Item Number: CM0003953



Introduction

KO validated Cyclin B1 rabbit monoclonal antibody for WB and ELISA applications, reacting with human. Targets a key cell cycle regulator essential for G2/M transition and CDK1 activation. Recombinant immunogen.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Cyclin B1 Rabbit mAb
Remarks/Alias	CCNB; B1
Species	Human
GeneID (human)	891
GeneID	891
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 3-177 of human Cyclin B1 (P14635)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P14635
Protein Weight	48kDa
Shipping	Ice bag