


[KO Validated] Anti-Cyclin D1 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P24385

Item Number: CM0003775



Introduction

KO-validated rabbit monoclonal antibody targeting Cyclin D1 (BCL1/PRAD1). Suitable for WB, IF/ICC, ELISA in human, mouse, and rat. Ideal for G1/S cell cycle transition and cancer research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Cyclin D1 Rabbit mAb
Remarks/Alias	BCL1; PRAD1; U21B31; D11S287E; D1
Species	Human
Gene ID (Human)	595
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 196-295 of human Cyclin D1 (P24385).
Source	Rabbit
Category	Monoclonal Antibody
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P24385
Protein Weight	34kDa
Shipping	Ice bag