

## [Ko Validated] Anti-Cytokeratin 20 (Krt20) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P35900

Item Number: CM0024403

### Introduction

High-quality rabbit polyclonal antibody targeting human cytokeratin 20 (CK20/KRT20), validated for WB, IF/ICC, and ELISA with cross-reactivity to mouse and rat. Ideal for studies on intestinal epithelial differentiation, keratin filament organization, and mucin secretion.

[Learn More](#)

| Parameter        | Value                                            |
|------------------|--------------------------------------------------|
| Product Name     | [KO Validated] Cytokeratin 20 (KRT20) Rabbit pAb |
| Remarks/Alias    | K20; CD20; CK20; CK-20; KRT21; 0)                |
| Species          | Human                                            |
| GeneID (human)   | 54474                                            |
| GeneID           | 54474                                            |
| Immunogen        | Recombinant protein                              |
| Source           | Rabbit                                           |
| Category         | Polyclonal Antibodies                            |
| Application      | WB, IF/ICC, ELISA                                |
| Cross Reactivity | Human, Mouse, Rat                                |
| SWISS            | P35900                                           |
| Protein Weight   | 48kDa                                            |
| Shipping         | Ice bag                                          |