


[KO Validated] Anti-Cytokeratin 8 Rabbit Polyclonal Antibody For Wb - P05787

Item Number: CM0019049



Introduction

KO validated rabbit polyclonal antibody targeting Cytokeratin 8 (KRT8). For Western blot applications, cross-reacts with human and mouse samples. Ideal for studying epithelial intermediate filaments, keratin networks, and related pathologies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Cytokeratin 8 Rabbit pAb
Remarks/Alias	K8; KO; CK8; CK-8; CYK8; K2C8; CARD2; 8
Species	Human
GeneID (human)	3856
Gene ID	3856
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 380-480 of human Cytokeratin 8 (NP_002264.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB
Cross Reactivity	Human, Mouse
SWISS	P05787
Protein Weight	54 kDa/56 kDa
Shipping	Ice bag