


Ko Validated Anti-Cytokeratin 8/Krt8 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P05787

Item Number: CM0009462



Introduction

KO-validated rabbit monoclonal antibody against human Cytokeratin 8/KRT8. Suitable for WB, IHC-P, and ELISA. Cross-reacts with human and rat. For epithelial biology and cancer research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Cytokeratin 8/KRT8 Rabbit mAb
Remarks/Alias	K8; CK8; CK-8; CYK8; K2C8; CARD2
Species	Human
GeneID (Human)	3856
GeneID	3856
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 341-440 of human Cytokeratin 8/KRT8 (NP_002264.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	P05787
Protein Weight	54 kDa/56 kDa
Shipping	Ice bag