

Ko Validated Anti-Ddit4 Polyclonal Antibody For Wb, Elisa - Q9Nx09

Item Number: CM0029641



Introduction

Rabbit polyclonal antibody targeting human DDIT4 (REDD1), KO validated for WB and ELISA. Ideal for studying mTOR signaling, cellular stress responses, and apoptosis in cancer and neurobiology research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] DDIT4 Rabbit pAb
Remarks/Alias	Dig2; REDD1; REDD-1; T4
Species	Human
GeneID	54541
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-232 of human DDIT4 (NP_061931.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9NX09
Protein Weight	25kDa
Shipping	Ice bag