


[Ko Validated] Anti-DJ-1/Park7 Polyclonal Antibody For Wb, If/Icc, Elisa - Q99497

Item Number: CM0013146



Introduction

Anti-DJ-1/PARK7 rabbit polyclonal antibody, validated for WB, IF/ICC, and ELISA. Reacts with human and mouse DJ-1/PARK7, a multifunctional oxidative stress sensor and deglycase. Ideal for cell protection and neurodegeneration research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] DJ-1/PARK7 Rabbit pAb
Remarks/Alias	DJ1; DJ-1; GATD2; HEL-S-67p; K7
Species	Human
GeneID (Human)	11315
GeneID	11315
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-189 of human DJ-1/PARK7 (NP_009193.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q99497
Protein Weight	20kDa
Shipping	Ice bag