

Ko Validated Anti-Dnajb1 Polyclonal Antibody For Wb, If/Icc, Elisa - P25685

Item Number: CM0025336



Introduction

KO-validated rabbit polyclonal antibody against DNAJB1 (Hdj1/HSP40). Suitable for WB, IF/ICC, and ELISA; cross-reacts with human, mouse, and rat. Ideal for chaperone and stress response studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] DNAJB1 Rabbit pAb
Remarks/Alias	Hdj1; Sis1; HSPF1; Hsp40; RSPH16B; B1
Species	Human
Gene ID (Human)	3337
Gene ID	3337
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-340 of human DNAJB1 (NP_006136.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P25685
Protein Weight	38kDa
Shipping	Ice bag