

[KO Validated] Anti-Ebp1/Pa2G4 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Uq80

Item Number: CM0016395



Introduction

Rabbit polyclonal antibody targeting EBP1/PA2G4, validated for WB, IHC-P, ELISA. Reacts with human, mouse, and rat. Suitable for research on cell proliferation, differentiation, and ERBB3 signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] EBP1/PA2G4 Rabbit pAb
Remarks/Alias	EBP1; HG4-1; ITAF45; p38-2G4; G4
Species	Human
Gene ID (Human)	5036
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 1-394 of human EBP1/PA2G4 (NP_006182.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UQ80
Protein Weight	44kDa
Shipping	Ice bag