

[KO Validated] Anti-Ece1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P42892

Item Number: CM0025546



Introduction

KO validated rabbit polyclonal antibody against human ECE1, suitable for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat ECE1. Ideal for endothelin-1 pathway research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] ECE1 Rabbit pAb
Remarks/Alias	ECE; E1
Species	Human
GeneID (human)	1889
GeneID	1889
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 670-770 of human ECE1 (NP_001388.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P42892
Protein Weight	87kDa
Shipping	Ice bag