

[KO Validated] Anti-Eef1B2 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P24534

Item Number: CM0029143



Introduction

CamelBio supplies a KO-validated rabbit polyclonal antibody targeting human EEF1B2, suitable for WB, IHC-P, IF/ICC, and ELISA in human, mouse, and rat. EEF1B2 functions as an elongation factor in protein biosynthesis.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] EEF1B2 Rabbit pAb
Remarks/Alias	EF1B; EEF1B; EEF1B1; B2
Species	Human
GeneID (Human)	1933
GeneID	1933
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-225 of human EEF1B2 (NP_001032752.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P24534
Protein Weight	25kDa
Shipping	Ice bag