

## Ko Validated Anti-Ehmt2 Polyclonal Antibody For Wb, Ip, Elisa - Q96Kq7

Item Number: CM0014248



### Introduction

KO validated rabbit polyclonal antibody against human EHMT2 (G9a), a histone H3K9 methyltransferase. Designed for Western blot (WB), immunoprecipitation (IP), and ELISA. Ideal for epigenetic and chromatin studies.

[Learn More](#)

| Parameter        | Value                                                                                                                                     |
|------------------|-------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | [KO Validated] EHMT2 Rabbit pAb                                                                                                           |
| Remarks/Alias    | G9A; BAT8; GAT8; NG36; KMT1C; C6orf30; T2                                                                                                 |
| Species          | Human                                                                                                                                     |
| Gene ID (Human)  | 10919                                                                                                                                     |
| Gene ID          | 10919                                                                                                                                     |
| Immunogen        | Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 877-1176 of human EHMT2 (NP_079532.5). |
| Source           | Rabbit                                                                                                                                    |
| Category         | Polyclonal Antibodies                                                                                                                     |
| Application      | WB, IP, ELISA                                                                                                                             |
| Cross Reactivity | Human                                                                                                                                     |
| SWISS            | Q96KQ7                                                                                                                                    |
| Protein Weight   | 132kDa                                                                                                                                    |
| Shipping         | Ice bag                                                                                                                                   |