


[KO Validated] Anti-Fascin/Fscn1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q16658

Item Number: CM0021166



Introduction

KO validated rabbit polyclonal antibody against human Fascin (FSCN1) for WB, IHC-P, IF/ICC, ELISA. Detects human, mouse, rat Fascin; suitable for actin cytoskeleton and cell motility research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Fascin/FSCN1 Rabbit pAb
Remarks/Alias	HSN; SNL; p55; FAN1; N1
Species	Human
Gene ID (Human)	6624
Gene ID	6624
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 394-493 of human Fascin/FSCN1 (NP_003079.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16658
Protein Weight	55kDa
Shipping	Ice bag