


[KO Validated] Anti-Galectin-3/Lgals3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P17931

Item Number: CM0012853



Introduction

KO validated rabbit polyclonal antibody against human Galectin-3 (LGALS3). Suitable for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Recognizes a 26kDa protein. Ideal for inflammation and cancer research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Galectin 3/LGALS3 Rabbit pAb
Remarks/Alias	L31; GAL3; MAC2; CBP35; GALBP; GALIG; LGALS2; S3
Species	Human
Gene ID (Human)	3958
Gene ID	3958
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human Galectin 3/LGALS3 (NP_002297.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17931
Protein Weight	26kDa
Shipping	Ice bag