

Ko Validated Anti-Gars Monoclonal Antibody For Wb, Elisa - P41250

Item Number: CM0003248



Introduction

KO validated rabbit monoclonal antibody against human GARS (glycyl-tRNA synthetase). Suitable for Western blot and ELISA. Immunogen corresponds to amino acids 640-739 of human GARS (P41250).

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] GARS Rabbit mAb
Remarks/Alias	GARS; HMN5; CMT2D; DSMAV; GlyRS; HMN5A; SMAD1; SMAJI; RS
Species	Human
Gene ID (Human)	2617
Gene ID	2617
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P41250
Protein Weight	83kDa
Shipping	Ice bag