

[KO Validated] Anti-Gata3 Rabbit Polyclonal Antibody For Wb And Elisa - P23771

Item Number: CM0019880



Introduction

Rabbit polyclonal antibody against human GATA3, validated by knockout. Suitable for WB and ELISA; cross-reacts with mouse and rat. For studying T-cell transcription factor function in Th2 differentiation and immune response.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] GATA3 Rabbit pAb
Remarks/Alias	HDR; HDRS; A3
Species	Human
Gene ID (Human)	2625
Gene ID	2625
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human GATA3 (NP_001002295.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23771
Protein Weight	48kDa
Shipping	Ice bag