

Ko Validated Anti-Git1 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Y2X7

Item Number: CM0020197

Introduction

KO validated rabbit polyclonal antibody against GIT1, a scaffold protein involved in focal adhesion, cell migration, and synaptic function. Suitable for WB, IF/ICC, ELISA in human, mouse, rat samples. Research use only.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] GIT1 Rabbit pAb
Remarks/Alias	p95-APP1; T1
Species	Human
Gene ID (Human)	28964
Gene ID	28964
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 460-640 of human GIT1 (NP_054749.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y2X7
Protein Weight	84kDa
Shipping	Ice bag