

[KO Validated] Anti-Glod4 Polyclonal Antibody For Wb, Elisa - Q9Hc38

Item Number: CM0022118



Introduction

High-quality rabbit polyclonal antibody against human GLOD4, validated for Western blot and ELISA. Suitable for human and mouse samples. Ideal for mitochondrial glyoxalase domain-containing protein 4 research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] GLOD4 Rabbit pAb
Remarks/Alias	HC6; HC71; CGI-150; C17orf25; D4
Species	Human
Gene ID (Human)	51031
Gene ID	51031
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 179-298 of human GLOD4 (NP_057164.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9HC38
Protein Weight	35kDa
Shipping	Ice bag