


KO Validated Anti-Glutathione Synthetase (Gss) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P48637

Item Number: CM0018508



Introduction

KO validated rabbit polyclonal antibody against human Glutathione Synthetase (GSS). Suitable for Western blot, IF/ICC, and ELISA; cross-reacts with mouse and rat. Ideal for oxidative stress studies and glutathione biosynthesis research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Glutathione Synthetase (GSS) Rabbit pAb
Remarks/Alias	GSHS; HEL-S-64p; HEL-S-88n; S)
Species	Human
GeneID (Human)	2937
GeneID	2937
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 373-474 of human Glutathione Synthetase (GSS) (NP_000169.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P48637
Protein Weight	52kDa
Shipping	Ice bag