

[KO Validated] Anti-Gsdme (Full Length+N Terminal) Monoclonal Antibody For Wb, Ip, Ihc-P, Elisa - O60443

Item Number: CM0009711



Introduction

High-quality [KO Validated] GSDME rabbit monoclonal antibody validated for WB, IP, IHC-P, and ELISA. Recognizes human GSDME (DFNA5/ICERE-1) with isoforms at 11, 35, and 55 kDa. Ideal for pyroptosis and apoptosis studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] GSDME (Full Length+N terminal) Rabbit mAb
Remarks/Alias	DFNA5; ICERE-1
Species	Human
GeneID (human)	1687
GeneID	1687
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 245-344 of human GSDME (Full Length+N terminal) (NP_004394.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IHC-P, ELISA
Cross Reactivity	Human
SWISS	O60443
Protein Weight	11 kDa/35 kDa/55 kDa
Shipping	Ice bag