

[KO Validated] Anti-Hadha Rabbit Polymab® For Wb, Ihc-P, If/Icc, Ip, Elisa - P40939

Item Number: CM0006414



Introduction

High-quality [KO Validated] anti-HADHA rabbit monoclonal antibody for WB, IHC-P, IF/ICC, IP, ELISA. Cross-reacts with human, mouse, rat. Essential for mitochondrial beta-oxidation and fatty acid metabolism research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] HADHA Rabbit PolymAb®
Remarks/Alias	GBP; ECHA; HADH; LCEH; MTPA; LCHAD; TP-ALPHA
Species	Human
GeneID (Human)	3030
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 545-763 of human HADHA (NP_000173.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P40939
Protein Weight	83 kDa
Shipping	Ice bag