

[KO Validated] Anti-Hdac1 Rabbit Monoclonal Antibody For Wb, Ip, If, Ihc, Chip, Elisa - Q13547

Item Number: CM0005700



Introduction

Rabbit monoclonal anti-HDAC1 antibody validated for WB, IP, IF, IHC, ChIP, ELISA. Suitable for human and mouse HDAC1 research in epigenetics and gene regulation.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] HDAC1 Rabbit mAb
Remarks/Alias	HD1; RPD3; KDAC1; GON-10; RPD3L1
Species	Human
GeneID (Human)	3065
Immunogen	Synthetic peptide (aa 393-482 of human HDAC1, NP_004955.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, IF-P, IHC-P, ChIP, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q13547
Protein Weight	55 kDa
Shipping	Ice bag