


[KO Validated] Anti-Her2/ErbB2 Polyclonal Antibody For Wb, Elisa - P04626

Item Number: CM0019307



Introduction

High-quality [KO Validated] anti-HER2/ErbB2 rabbit polyclonal antibody validated for Western blot (WB) and ELISA. Targets human HER2/ErbB2, a receptor tyrosine kinase overexpressed in many cancers. Ideal for oncology research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] HER2/ErbB2 Rabbit pAb
Remarks/Alias	NEU; NGL; HER2; TKR1; CD340; HER-2; VSCN2; MLN 19; MLN-19; c-ERB2; c-ERB-2; HER-2/neu; p185(erbB2); B2
Species	Human
GeneID	2064
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1156-1255 of human HER2/ErbB2 (NP_004439.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P04626
Protein Weight	138kDa
Shipping	Ice bag