

[KO Validated] Anti-Hexokinase-1 (Hk1) Rabbit Monoclonal Antibody For Wb, Ip, Ihc-P, Elisa - P19367

Item Number: CM0002664



Introduction

KO-validated rabbit monoclonal antibody targeting human Hexokinase-1 (HK1, P19367). Suitable for WB, IP, IF/ICC, IF-P, IHC-P, ELISA; reacts with human, mouse, rat. Supports glycolysis research, metabolism, and innate immunity studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] HK1 Rabbit mAb
Remarks/Alias	HK; HKD; HKI; HXK1; NMSR; RP79; HMSNR; HK1-ta; HK1-tb; HK1-tc; NEDVIBA; hexokinase; K1
Species	Human
GeneID (Human)	3098
GeneID	3098
Immunogen	Synthetic Peptide (A synthetic peptide corresponding to a sequence within amino acids 50-150 of human Hexokinase 1 (P19367).)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P19367
Protein Weight	102kDa
Shipping	Ice bag