

Ko Validated Anti-Hexokinase Ii Monoclonal Antibody For Wb, Elisa - P52789

Item Number: CM0007708



Introduction

KO validated Rabbit monoclonal anti-Hexokinase II antibody for Western blot and ELISA; detects human, mouse, and rat Hexokinase II (102 kDa). Ideal for glycolysis and mitochondrial apoptosis studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Hexokinase II Rabbit mAb
Remarks/Alias	HKII; HXK2; II
Species	Human
GeneID (human)	3099
GeneID	3099
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human Hexokinase II (NP_000180.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P52789
Protein Weight	102kDa
Shipping	Ice bag