


[KO Validated] Anti-Hira Polyclonal Antibody For Wb, Ihc-P, Elisa - P54198

Item Number: CM0030591



Introduction

Rabbit polyclonal antibody against human HIRA (TUP1/DGCR1), KO validated for WB, IHC-P, and ELISA applications. Useful for studying chromatin assembly, senescence, and gene regulation.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] HIRA Rabbit pAb
Remarks/Alias	TUP1; DGCR1; TUPLE1; RA
Species	Human
GeneID (Human)	7290
GeneID	7290
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 250-550 of human HIRA (NP_003316.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P54198
Protein Weight	112kDa
Shipping	Ice bag