


[KO Validated] Anti-Histone H2A.z Polyclonal Antibody For Wb, Elisa - P0C0S5

Item Number: CM0029776



Introduction

Rabbit polyclonal antibody against human Histone H2A.Z, KO validated for WB and ELISA. Cross-reacts with mouse and rat. Immunogen: recombinant protein (aa 1-128). Molecular weight ~13 kDa. Ideal for chromatin biology and epigenetic research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Histone H2A.Z Rabbit pAb
Remarks/Alias	H2AZ; H2A.z; H2A/z; H2AFZ; H2A.Z-1; .Z
Species	Human
Gene ID (Human)	H2A
Gene ID	3015
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-128 of human H2AFZ (NP_002097.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P0C0S5
Protein Weight	13 kDa
Shipping	Ice bag