


Ko Validated Anti-Hmgb1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P09429

Item Number: CM0029892



Introduction

CamelBio offers KO validated anti-HMGB1 rabbit polyclonal antibody (P09429) for Western blot, IHC-P, and ELISA. Recognizes human, mouse, rat HMGB1. Ideal for DAMP and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] HMGB1 Rabbit pAb
Remarks/Alias	HMG1; HMG3; HMG-1; SBP-1; HMGB1
Species	Human
GeneID (Human)	3146
GeneID	3146
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 69-168 of human HMGB1 (NP_002119.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09429
Protein Weight	25 kDa
Shipping	Ice bag