

[KO Validated] Anti-Hnrnp E1/Pcbp1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q15365

Item Number: CM0018959



Introduction

Rabbit polyclonal antibody targeting human hnRNP E1/PCBP1, validated for Western blot, IHC-P, and ELISA. Suitable for human, mouse, and rat samples, with a molecular weight of 37 kDa. Ideal for studying RNA-binding protein functions, erythropoiesis, and viral infection mechanisms.

[Learn More](#)

Parameter	Specification
Product Name	[KO Validated] hnRNP E1/PCBP1 Rabbit pAb
Remarks/Alias	HNRPX; HNRPE1; hnRNP-X; HEL-S-85; hnRNP-E1; P1
Species	Human
Gene ID	5093
Immunogen	Synthetic peptide; A synthetic peptide corresponding to a sequence within amino acids 100-200 of human hnRNP E1/PCBP1 (NP_006187.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15365
Protein Weight	37kDa
Shipping	Ice bag