

[KO Validated] Anti-Hnrnpa2B1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P22626

Item Number: CM0013846



Introduction

[KO Validated] Rabbit polyclonal antibody against HNRNPA2B1 for WB, IHC-P, IF/ICC, ELISA. Detects human, mouse, rat. Ideal for RNA biology and telomere research. - P22626

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] HNRNPA2B1 Rabbit pAb
Remarks/Alias	RNPA2; HNRPA2; HNRPB1; SNRPB1; HNRNPA2; HNRNPB1; IBMPFD2; HNRPA2B1; B1
Species	Human
Gene ID (Human)	3181
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human HNRNPA2B1 (NP_112533.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P22626
Protein Weight	37kDa
Shipping	Ice bag