


[KO Validated] Anti-Hprt1 Rabbit Polyclonal Antibody For Wb, Elisa - P00492

Item Number: CM0028117



Introduction

Rabbit polyclonal antibody targeting HPRT1, validated for knockout applications. Suitable for WB and ELISA in human, mouse, and rat samples, supporting purine salvage pathway research.

[Learn More](#)

Parameter	Details
Product Name	[KO Validated] HPRT1 Rabbit pAb
Remarks/Alias	HPRT; HGPRT; T1
Species	Human
GeneID (Human)	3251
GeneID	3251
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-218 of human HPRT1 (NP_000185.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00492
Protein Weight	25kDa
Shipping	Ice bag