

[KO Validated] Anti-Hsp90A Monoclonal Antibody For Wb, Elisa - P07900

Item Number: CM0011650



Introduction

Rabbit monoclonal antibody against human Hsp90 α , validated for WB and ELISA. Cross-reacts with mouse and rat. Suitable for chaperone biology, signal transduction, and inflammatory response studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Hsp90 α Rabbit mAb
Remarks/Alias	EL52; HSPN; LAP2; HSP86; HSPC1; HSPCA; Hsp89; Hsp90; LAP-2; HSP89A; HSP90A; HSP90N; Hsp103; HSPCAL1; HSPCAL4; HEL-S-65p
Species	Human
Gene ID	3320
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 433-732 of human Hsp90 α (NP_005339.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07900
Protein Weight	84kDa/98kDa
Shipping	Ice bag