


[Ko Validated] Anti-Hsp90B Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P08238

Item Number: CM0019332



Introduction

Rabbit polyclonal antibody against human Hsp90 β (HSP90AB1), validated for KO, suitable for WB, IF/ICC, ELISA applications. Reacts with human, mouse, rat, monkey. Molecular chaperone involved in protein folding and signal transduction.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Hsp90 β Rabbit pAb
Remarks/Alias	HSP84; HSPC2; HSPCB; D6S182; HSP90B; 0 β
Species	Human
GeneID (Human)	3326
GeneID	3326
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, Monkey
SWISS	P08238
Protein Weight	83kDa
Shipping	Ice bag