

[KO Validated] Anti-Idh1 Rabbit Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - O75874

Item Number: CM0013137



Introduction

[KO Validated] rabbit polyclonal antibody against human IDH1, suitable for WB, IF/ICC, IP, and ELISA applications. Detects IDH1 in human, mouse, and rat samples, targeting the cytoplasmic NADP-isocitrate dehydrogenase.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] IDH1 Rabbit pAb
Remarks/Alias	IDH; IDP; IDCD; IDPC; PICD; HEL-216; HEL-S-26; H1
Species	Human
GeneID (Human)	3417
GeneID	3417
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-414 of human IDH1 (NP_005887.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75874
Protein Weight	47kDa
Shipping	Ice bag