

[KO Validated] Anti-Igf2Bp1/Imp1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Nzi8

Item Number: CM0019809



Introduction

KO validated rabbit polyclonal antibody against human IGF2BP1/IMP1. Suitable for WB and ELISA, cross-reacts with mouse and rat. Targets an RNA-binding protein involved in mRNA transport, stability, and translation regulation in cancer and neurobiology.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] IGF2BP1/IMP1 Rabbit pAb
Remarks/Alias	IMP1; ZBP1; CRDBP; IMP-1; CRD-BP; VICKZ1; P1
Species	Human
Gene ID (Human)	10642
Immunogen	Recombinant fusion protein containing amino acids 1-100 of human IMP1 (NP_006537.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NZI8
Protein Weight	63kDa
Shipping	Ice bag