


[KO Validated] Anti-Impdh2 Monoclonal Antibody For Wb, If/Icc, Elisa - P12268

Item Number: CM0008005



Introduction

KO-validated rabbit monoclonal antibody against human IMPDH2, suitable for WB, IF/ICC, and ELISA. Detects human, mouse, and rat IMPDH2; key for nucleotide metabolism research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] IMPDH2 Rabbit mAb
Remarks/Alias	IMPD2; IMPDH-II; H2
Species	Human
GeneID	3615
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-300 of human IMPDH2 (P12268)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P12268
Protein Weight	56 kDa
Shipping	Ice bag