


[KO Validated] Anti-Ing5 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q8WYh8

Item Number: CM0029343



Introduction

Rabbit polyclonal antibody against human ING5, validated for Western blot, IF/ICC, and ELISA. Suitable for chromatin acetylation, cell cycle, and apoptosis research. Cross-reacts with human.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] ING5 Rabbit pAb
Remarks/Alias	p28ING5; G5
Species	Human
GeneID (Human)	84289
GeneID	84289
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human ING5 (NP_115705.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q8WYH8
Protein Weight	28kDa
Shipping	Ice bag