

[KO Validated] Anti-Ipo4 Rabbit Polyclonal Antibody For Wb, Elisa - Q8Tex9

Item Number: CM0019615



Introduction

KO-validated rabbit polyclonal antibody targeting human Importin-4 (IPO4) for WB and ELISA. Cross-reacts with mouse and rat. Key reagent for investigating nuclear transport of histones, VDR, and viral proteins.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] IPO4 Rabbit pAb
Remarks/Alias	Imp4; O4
Species	Human
GeneID	79711 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 728-1081 of human IPO4 (NP_078934.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8TEX9
Protein Weight	119kDa
Shipping	Ice bag