

[KO Validated] Anti-Irf3 Rabbit Monoclonal Antibody For Wb - Q14653

Item Number: CM0002222



Introduction

KO validated rabbit monoclonal antibody targeting interferon regulatory factor 3 (IRF3). Suitable for Western blot (WB) and ELISA applications. Cross-reacts with human and mouse IRF3 protein. Ideal for innate immune response and antiviral research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] IRF3 Rabbit mAb
Remarks/Alias	IIAE7; F3
Species	Human
Gene ID	3661
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 50-150 of human IRF3 (Q14653).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q14653
Protein Weight	47kDa
Shipping	Ice bag