

[Ko Validated] Anti-Irf9 Polyclonal Antibody For Wb, Elisa - Q00978

Item Number: CM0021460

Introduction

KO-validated rabbit polyclonal antibody targeting human IRF9 (ISGF3γ), suitable for WB and ELISA. Detect the 44 kDa transcription factor crucial for type I interferon antiviral signaling. For innate immunity and ISGF3 complex studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] IRF9 Rabbit pAb
Remarks/Alias	p48; IRF-9; ISGF3; ISGF3G; F9
Species	Human
Gene ID	10379
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 113-255 of human IRF9 (NP_006075.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q00978
Protein Weight	44kDa
Shipping	Ice bag