


Ko Validated Anti-Itch Polyclonal Antibody For Wb, Elisa - Q96J02

Item Number: CM0016045



Introduction

KO Validated rabbit polyclonal antibody against ITCH (E3 ubiquitin ligase Itchy homolog) for WB, ELISA in human, mouse, rat. Involved in ubiquitination, inflammation, immune response, and apoptosis.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] ITCH Rabbit pAb
Remarks/Alias	AIF4; AIP4; ADMFD; NAPP1; CH
Species	Human
Gene ID (Human)	83737
Gene ID	83737
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-155 of ITCH (NP_001244066.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96J02
Protein Weight	103kDa
Shipping	Ice bag